





















### INDICATIONS

Nevus (birth mark)
Seborrheic keratosis
Wart
Corn, Clavus
Acne
Circumcision
Operation on snoring
Hemorrhoid
Lentigo
Milium
Syringoma
Soft fibroma
Ingrowing nail
Wrinkle removal
Axillary odor, etc.

# 15 INDICATION MODES





## 15 Indicaion modes

In the "indication mode", the energy and procedural guidance for each symptom is indicated on the display, offering practitioners a quick and easy procedural environment. In addition, the practitioners can control "repeat time", "pulse width", etc., according to the condition of the patient, enabling an even more detailed and stable procedure.

#### General

Dimensions (w × h × d) Weight Weight including accessories

Main voltage Maximum power, in operation

Safety class Fuses

Aimming beam

Cooling system

320 × 1160 × 420 mm

40.3 kg 44.1 kg

220~240 VAC, 50~60 Hz

300 VA

Class IV laser according to IEC 60825-1

2 × Little AC 250V 10A

#### Laser module

Wavelength
Beam mode quality
Beam quality (at 1/e²)
Laser power output, continuous
Dot size
Delivery system

10.57 ~ 10.63  $\mu$ m, invisible infrared TEMoo, 95% purity

M<sup>2</sup>≤ 1.2

30 Watts at arm's tail

7-joint, fixed mirror articulated arm 5mW, 650nm red diode laser

Water

#### UltraPulse mode

Pulse width  $200 \sim 1000 \, \mu \text{s}$ Repeat time  $10 \sim 1000 \, \text{ms}$ 

#### SuperPulse mode

Pulse width  $1 \sim 20 \text{ ms}$ Repeat time  $50 \sim 1000 \text{ ms}$ 

#### Continuous Wave mode

Power 0.5 ~ 30 W

#### **Environmental Condition**

Temperature  $+10 \sim +40 \,^{\circ}\text{C}$ Relative humidity  $30 \sim 75\%$ Atmospheric pressure  $700 \sim 1060 \text{ hPa}$ 



BISON Medical CO.,LTD.

www.bisonmedical.com

205 Seungil Venture Tower, 371-42 Gasan-dong, Geumcheon-gu, Seoul 153-803, Korea
Tel) +82-2-865-7121
Fax) +82-2-865-7131
sales@bisonmedical.com