

DentalQ

Dental Diode Laser

Selectable Among
808 / 915 / 940 / 976 nm

- Selectable Wavelength
- Compact Size Device With Hand Grip
- Multifunctional Handpieces With Fiber
- Intuitive Handling
- Touch Screen Interface




DentalQ Indication

- | | |
|---|---|
|  1. Pulp capping |  8. Hypersensitivity _ occlusal and cervical |
|  2. Tooth canal irrigation |  9. Dental caries |
|  3. Periodontal abscess |  10. Gingival incision / Gingivoplasty |
|  4. Hemostasis after tooth extraction |  11. Periodontal pocket irrigation |
|  5. Fluoride application & Tooth bleaching |  12. Periodontal disease |
|  6. Removal of tooth pigmentation |  13. User Indication |
|  7. Implant surgery | |

12 Indications + 1 User Custom

Save 30 on Indication 1 = Total 390 slots



General

Dimensions (W x D x H)	167(w) x 146.3(D) x 236(H) (mm)
Weight	1.6 kg
Main voltage	12VDC Adaptor
Safety class	Class IV laser
Connection	SMA905

Laser module

Wavelength	Selectable Among 808 / 915 / 940 / 976 nm
Laser power output, continuous	808nm / 7W 915nm / 8W 940nm / 8W 976nm / 8W
Delivery system	Fiber patch cord
Fiber core diameter	200 / 300 / 400μm
Aiming beam	Red 650nm, ~5mW

B I S O N
est
deal
olution of
optical
ano-technology

Seungil Venture Tower, 136, Gasan digital 2-ro, Geumcheon-gu, Seoul, Korea
Office: 503(5F), R&D: 504(5F), Factory: 205(2F), 804(8F)

Tel) +82-2-865-7121 Fax) +82-2-865-7131

dental@bisonmedical.com

Copyright© BISON MEDICAL CO., LTD. All rights reserved.

BISON
MEDICAL EQUIPMENT MANUFACTURER